

MEDTRONIC IS COMMITTED TO THE BRAIN MODULATION BUSINESS

We alleviate pain, restore health and extend life by delivering a Medtronic brain modulation solution to every eligible patient.

Off Label Disclosure Slide

Thank you for participating in this educational event!


This program, sponsored by Medtronic, is intended to *educate* and *train* customers regarding the **approved** or **cleared** uses of Medtronic products. As such, unapproved products or indications may not be presented or discussed during the program.

You may contact Medtronic's Office of Medical Affairs at 800.876.3133 ext. 6044 or rs.msdoma@medtronic.com for any specific clinical questions you may have.

MEDTRONIC DBS THERAPY— IMPROVING LIVES THEN, NOW, TOMORROW.

1980s 1990s 2000s 2010s

1987
Medtronic DBS Therapy is implanted for the first time in Grenoble, France.



1995
10,000th person worldwide implanted with Medtronic DBS Therapy.

CE Mark approval in Europe for treating essential tremor and tremor associated with Parkinson's disease.

1997
Medtronic DBS Therapy receives Food and Drug Administration (FDA) approval to treat essential tremor and tremor associated with Parkinson's disease.



1998
Medtronic DBS Therapy approved in Europe for Parkinson's disease.



2002
Medtronic receives FDA approval of DBS Therapy for the treatment of Advanced PD.

Advanced PD

2003
Medtronic DBS Therapy receives approval in Europe and Humanitarian Device Exemption (HDE)* approval in the U.S. for managing the symptoms of primary dystonia.



2007
50,000th person worldwide implanted with Medtronic DBS Therapy.



2009
Medtronic DBS Therapy receives approval in Europe and an HDE* approval in the U.S. for the treatment of Obsessive-Compulsive Disorder (OCD).



2010
European approval of Medtronic DBS Therapy to treat medically refractory epilepsy patients with partial-onset seizures.†

2012
Medtronic DBS Therapy receives approval in Australia and Canada for the treatment of medically refractory epilepsy patients with partial-onset seizures.†

2015
CE Mark in Europe and FDA approval of full-body MR Conditional MRI scans under specific conditions. FDA Approval of Expanded Labeling.



2016
140,000th person worldwide implanted with Medtronic DBS Therapy.

US Launch of Expanded PD Labeling

2017
160,000th person worldwide implanted with Medtronic DBS Therapy

2018
Medtronic DBS Therapy receives Food and Drug Administration (FDA) approval to treat as an **adjunctive treatment for reducing the frequency of partial-onset seizures in those who are refractory to ≥3 antiepileptic medications.**




*Humanitarian device exemption (HDE) in the United States; the effectiveness of this device for the treatment of dystonia and obsessive-compulsive disorder (OCD) has not been demonstrated.

DBS for Refractory Epilepsy



IMPLANTABLE SYSTEM COMPONENTS



Three Components*

1. Implantable Neurostimulator (INS): Power
2. Extension: connects the INS to the lead
3. Lead: Implanted in the brain, electrodes in contact with target tissue

*Some systems may include a pocket adaptor


MEDTRONIC DBS SYSTEM FOR EPILEPSY

- Activa PC Model 37601 Neurostimulator



- Model 3387/3389 DBS Lead

- The DBS Leads have four electrodes



- Model 37086 Extension

- Distal (Extension connection) end has four contacts
- Available in lengths
 - 40 cm
 - **60 cm**
 - 95 cm
- Stretch coil extension allows for up to 15% extensibility



INTERCEPT™ PATIENT PROGRAMMER

PC ONLY




HOW TO TURN NEUROSTIMULATOR ON AND OFF



MODE THERAPY SCREEN


▪ Simple Mode Therapy Screen shows:

- Therapy ON/OFF status
- Seizure Count



▪ Advanced Mode Therapy Screen shows:

- Therapy ON/OFF status
- Seizure Count
- Parameter settings and active group



INTERCEPT™ PATIENT PROGRAMMER MODES


- There are three modes available
 - Simple Mode
 - Advanced View
 - Advanced Adjust
- Mode is determined by clinician based on patient needs

	Simple Mode	Advanced View	Advanced Adjust
Turn therapy ON/OFF	•	•	•
Check INS Battery Status	•	•	•
Daily Battery Check Reminder	•	•	•
Additional Patient Programmer Features	•	•	•
View Therapy Parameters		•	•
Change Groups		•	•
Ability to adjust stimulation within the patient limits that are set for one or more group(s).			•

SETTING SEIZURE MODE

- Turn patient programmer ON (do not interrogate)
- Press and hold Selection Keys until Lead Connections Screen Appears
- Navigate to Seizure Mode Screen
- Set to **on** or **off**
- This is meant to be utilized by HCP programmers only and **NOT** the patient

Note: Default setting for Seizure mode is off



DBS THERAPY FOR EPILEPSY

INITIAL PROGRAMMING ACTIVITIES


- Configure leads (completed in SET-UP)
- Verify system integrity (check electrode impedances)
- Program initial stimulation parameters
- Program the neurostimulator for patient control
- Provide patient and caregiver with instructions on use of patient programmer and tracking of seizures (count, type, severity)
- Emphasize adherence to AED regimen
- Verify tolerability of stimulation

SANTE STIMULATION PARAMETERS: EPILEPSY

Parameter	Typical Starting Value
Amplitude	5 V *In clinical practice clinicians may start at a lower amplitude and titrate up slower as they assess patient response and tolerability.
Pulse Width	90 μ s
Rate	145 Hz
Electrode Configuration	Unipolar Mode: Single electrode or two adjacent electrodes negative, case positive (all patients in the SANTE clinical trial were in unipolar mode)
Cycle of Therapy	Cycling mode ON: 1 minute on, 5 minutes off
SoftStart™Stop	programmed to 8 seconds

Tablet Clinician Programmer

Launched June 2018



Modern, usable interface.

Focus on control of stimulation and documentation of outcomes.

Intuitive patient management.

A SMARTER DBS JOURNEY


PLATFORM OVERVIEW

A new era of innovation from the only partner with a **25-year DBS legacy**, **the Medtronic clinician programmer** offers easy navigation and an intuitive workflow over a safe and secure platform, so you can connect with patients with comfort and confidence.

STREAMLINING WORKFLOWS
Intuitive interface means more time focusing on your patient.


ENABLING SECURE PATIENT CONNECTIONS
Freedom to move with peace of mind.

INFORMATION MADE ACTIONABLE
Annotations and historical programming information at a glance that are saved in the patient's INS.



CLINICIAN PROGRAMMER & ACCESSORIES

COMMUNICATOR AND DRAPE



- * Optional protective case and drape may be ordered from Medtronic
- ** Activa App is NOT compatible with Soletra/Kinetra Neurostimulators

SYSTEM COMPONENTS

1.) TABLET *(with Activa App)

- Compatible with ALL Activa™ Family devices**

2.) COMMUNICATOR

- Encrypted Bluetooth connection from the programmer to the communicator
- Proprietary, proximal telemetry from the communicator to the implanted device

3.) DRAPE*

- Freedom of movement


INTENDED USE

- THIS DECK IS INTENDED TO BE USED TO TRAIN CLINICIANS ON THE MEDTRONIC CLINICIAN PROGRAMMER**
- IT IS NOT INTENDED TO BE USED A SUBSTITUTE OR REPLACEMENT FOR THE LABELING**
- INFORMATION CONTAINED HEREIN IS THE PROPERTY OF MEDTRONIC**
- ALL RIGHTS RESERVED ALL RIGHTS RESERVED. MEDTRONIC, MEDTRONIC LOGO AND FURTHER, TOGETHER ARE TRADEMARKS OF MEDTRONIC. ALL OTHER BRANDS ARE TRADEMARKS OF A MEDTRONIC COMPANY**

Medtronic
Further, Together

STREAMLINED

MODERN TECHNOLOGY
UNLEASHED



STREAMLINED

MODERN TECHNOLOGY
UNLEASHED

DEMO: CONNECT

SELECT DEVICE...


- Activa RC – new
- Activa RC
- Activa PC – new
- Activa PC
- Activa SC – new



CANCEL

RETRY



CONNECT



STREAMLINED

MODERN TECHNOLOGY
UNLEASHED

- ~:0/5 - CYCLING ON
- ~X:5/5 - CYCLING OFF



PATIENT

Patient Name: [No name]
Patient ID:
Diagnosis:

DEVICE

Device Model: Activa PC
Model Number: 37601
Serial Number: NKM728744
Implant Date: Oct 18, 2018
Battery Level: 2.95 V | OK

IMPEDANCE

Status: Perform an electrode impedance measurement

CLINICIAN NOTES

~:0/5 (0) ← ~ = Cycling ON
~X:5/5 (2) ← ~X = Cycling OFF

SELECT TASK

- SETUP
- STIMULATION
- IMPEDANCE
- REPLACEMENT
- END SESSION


What this means:


The patient hit the BLUE button on their programmer 5 times when the SEIZURE MODE was off (just an EVENT counter)

2 of the (5 times) times they hit the BUTTON the DBS therapy was COMPLETELY OFF.






NUMBER in Brackets **(2)**


STREAMLINED TASK BASED WORKFLOWS




Stimulation On 

SELECT TASK

- SETUP** 
- STIMULATION 
- IMPEDANCE 
- REPLACEMENT 
- END SESSION 






DEMO: HOME | John Smith | Stimulation On 

PATIENT
Patient Name: John Smith
Patient ID: MRN:123456
Diagnosis: Epilepsy

DEVICE
Device Model: Activa PC
Model Number: 37601
Serial Number: NKM_FOLLOWUP_PC
Implant Date: Apr 22, 2017
Battery Level:  3.81 V | OK
Cycling: Enabled
On Duration: 1 Minute
Off Duration: 5 Minutes

IMPEDANCE
Status: INVESTIGATE
LEFT HEMISPHERE


SELECT TASK

- SETUP 
- STIMULATION 
- IMPEDANCE 
- REPLACEMENT 
- END SESSION 

NO ALERTS

STREAMLINED

**EASY TO
NAVIGATE**



Stimulation On

SELECT TASK

- SETUP
- STIMULATION**
- IMPEDANCE
- REPLACEMENT
- END SESSION

DEMO: HOME | John Smith | Stimulation On

SELECT TASK

- SETUP
- STIMULATION
- IMPEDANCE
- REPLACEMENT
- END SESSION

PATIENT
Patient Name: John Smith
Patient ID: MRN:123456
Diagnosis: Epilepsy


DEVICE
Device Model: Activa PC
Model Number: 37601
Serial Number: NKM_FOLLOWUP_PC
Implant Date: Apr 22, 2017
Battery Level: 3.81 V | OK
Cycling: Enabled
On Duration: 1 Minute
Off Duration: 5 Minutes

IMPEDANCE
Status: INVESTIGATE
LEFT HEMISPHERE

NO ALERTS

INSIGHTFUL INFORMATION YOU NEED

Group & Program 'Management



Left ANT

Right ANT

DRAG AND DROP GROUPS

GROUP A

5.0 V
90 μ s
145 Hz

+

2.7 V
60 μ s
145 Hz

+

GROUP B

2.0 V
90 μ s
145 Hz

0.0 V
90 μ s
145 Hz

GROUP C

5.0 V
90 μ s
125 Hz

1.1 V
90 μ s
125 Hz

2.7 V
60 μ s
125 Hz

GROUP D

5.0 V
90 μ s
145 Hz

5.0 V
90 μ s
145 Hz

Select a group to activate or a program to edit. Drag items to copy or delete.



NO ALERTS

VIEW HISTORY

VIEW HISTORY


GROUPS

- “Packages” of therapy settings (electrode configurations, parameters, cycles)
- Clinicians may assign up to four groups using the clinician programmer (A, B, C, D) from which patients can select using their patient programmer
- A Group can be set in bipolar for MRI if conditions are met



STREAMLINED

SEE
BRIGHTER



Left ANT **LEAD LOCATION** SETTINGS

Lead Annotation

NAVIGATION TABS

VNA

ELECTRODE CONFIGURATION

PROGRAM TABS

CONTROLS

AMPLITUDE RESOLUTION

PARAMETER TABS

STIMULATION

Program 1 Program 2

AMPLITUDE

5.0 Volts


5.0 V 90 μ s 145 Hz

PATIENT LIMITS

NO ALERTS

SETTINGS

ADVANCED PROGRAMMING FEATURES



Left ANT

Lead Annotation

3
2
1
0

NO ALERTS

SETTINGS

STIMULATION

Program 1 Program 2

AMPLITUDE


5.0
Volts

5.0 V 90 μ s 145 Hz

PATIENT LIMITS


STREAMLINED

APP BASED
CONVENIENCE



The screenshot shows the 'DEMO: SETTINGS' screen for a user named John Smith. At the top right, there is a 'Stimulation' toggle switch which is turned 'On'. The settings are organized into sections: 'Stimulation' (Ramp Interval: 2 Seconds, Amplitude Control: Voltage, SoftStart/Stop: 4 Seconds, Cycling: On Duration: 1 Minute, Off Duration: 5 Minutes), 'Patient Programmer' (Mode: Advanced Adjust, Adjustable Parameter: Amplitude), and 'Check Battery Reminder' (11:00 AM). Each of these sections has a corresponding green 'On' toggle switch. A 'CLOSE' button is located at the bottom right of the screen. The top status bar shows 84% battery and 9:37 AM.

STREAMLINED
APP BASED
CONVENIENCE



Stimulation

Stimulation

Ramp Interval
2 Seconds

Amplitude Control
Voltage

SoftStart/Stop
4 Seconds

Cycling
On Duration: 1 Minute Off Duration: 5 Minutes

Patient Programmer

Mode
Advanced Adjust

Adjustable Parameter
Amplitude



Check Battery Reminder
11:00 AM

CLOSE

STIMULATION CYCLE OPTIONS SETTING UP ON/OFF DURATION

- Minimum off or on time: 0.1sec
- Maximum off or on time: 24 hours
- SANTE Starting Parameters-Intervals:

1 min on 5 Minutes off





* Enabling cycling at certain parameter settings may decrease the device longevity of non-rechargeable devices.

STREAMLINED

APP BASED
CONVENIENCE

Setting Patient Limits



Stimulation

Ramp Interval
2 Seconds

Amplitude Control
Voltage

SoftStart/Stop
4 Seconds

Cycling
On Duration: 1 Minute Off Duration: 5 Minutes

Patient Programmer

Mode
Advanced Adjust

Adjustable Parameter
Amplitude

Check Battery Reminder
11:00 AM

CLOSE


LETS PRACTICE SETTING UP A SAMPLE- PC- WORK IN DEMO MODE
 NOT NEW – IF YOU HAVE ALL YOUR GROUPS FILLED DELETE THEM.

Parameter	Value
Amplitude	2.5 V
Pulse Width	90 μ s
Rate	145 Hz
Electrode Configuration	2 (-) negative, case (+)
Cycle of Therapy	Cycling mode ON: 1 minute on, 5 minutes off
SoftStart™Stop	programmed to 8 seconds
Set up NEW for MRI	2(-) 3 (+) BIPOLAR so patient may be able to stay on during the MRI

- Final task – What are the actions you want to take at the end of the session
- Be sure you end the session with the group you want the patient to leave the office on.

STREAMLINED

EASY TO
NAVIGATE



Stimulation On

SELECT TASK

- SETUP
- STIMULATION
- IMPEDANCE**
- REPLACEMENT
- END SESSION

Stimulation On

DEMO: HOME | John Smith

PATIENT
Patient Name: John Smith
Patient ID: MRN:123456
Diagnosis: Epilepsy

DEVICE
Device Model: Activa PC
Model Number: 37601
Serial Number: NKM_FOLLOWUP_PC
Implant Date: Apr 22, 2017
Battery Level: 3.81 V | OK
Cycling: Enabled
On Duration: 1 Minute
Off Duration: 5 Minutes

IMPEDANCE
Status: INVESTIGATE
LEFT HEMISPHERE

SELECT TASK

- SETUP
- STIMULATION
- IMPEDANCE
- REPLACEMENT
- END SESSION

NO ALERTS

**INSIGHTFUL
INFORMATION
MADE
ACTIONABLE**

You can select an automated impedance test or select a test amplitude for the electrode impedance measurement.

If the automated selection is chosen the neurostimulator will automatically use the following combinations of stimulation parameters to measure impedance:

	Measurement Output (V) amplitude			Measurement Range	Measurement At High Output
	PW	Rate			
Meas 1	80µs	100 Hz	0.7	Up to 10,000 ohms	Value above 2,000
Meas 2	80µs	100 Hz	1.5	Up to 20,000 ohms	Value above 2,000
Meas 3	80µs	100 Hz	3.0	Up to 40,000 ohms	No prompt: highest range already achieved

-The system takes the measurements in sequence. It is possible that the system will not need to take all three measurements.

DEMO: IMPEDANCE | John Smith | Stimulation On

Summary **Electrode** Therapy

Electrode impedance tests verify the electrical integrity of individual electrodes (monopolar) and pairs of electrodes (bipolar).

Left ANT

	MONOPOLAR	BIPOLAR		
		0	1	2
3	✓ 1,270	! 15.2K	✓ 1,590	✓ 2,262
2	✓ 1,310	! 15.7K	✓ 1,865	
1	✓ 1,319	! 15.4K		
0	! 14.0K			

GUIDANCE

- ✓ OK: Impedance is within normal range.
- ! INVESTIGATE: Impedance is outside normal range.
- ✗ AVOID: Impedance indicates a possible fault condition.

MEASURE ELECTRODE IMPEDANCE

ALERTS: 1

CLOSE

STREAMLINED

EASY TO
NAVIGATE

Stimulation



SELECT TASK

SETUP



STIMULATION



IMPEDANCE



REPLACEMENT



END SESSION



DEMO: HOME | John Smith

Stimulation On

PATIENT
Patient Name: John Smith
Patient ID: MRN:123456
Diagnosis: Epilepsy


DEVICE
Device Model: Activa PC
Model Number: 37601
Serial Number: NKM_FOLLOWUP_PC
Implant Date: Apr 22, 2017
Battery Level: 3.81 V | OK
Cycling: Enabled
On Duration: 1 Minute
Off Duration: 5 Minutes

IMPEDANCE
Status: INVESTIGATE
LEFT HEMISPHERE

SELECT TASK

- SETUP
- STIMULATION
- IMPEDANCE
- REPLACEMENT
- END SESSION

NO ALERTS



STREAMLINED

SECURE

INSIGHTFUL

REVIEWING FINAL SETTINGS AND ENDING A SESSION

Review Final Settings...

1

- Swipe up and down to scroll
- Left side = initial settings
- Right side = final settings

2

- Undo button = to undo programming changes made during programming session

3

Go to Stimulation = go directly to Stimulation screen and adjust stimulation settings further

4

Clinician Notes = enter free-text notes

5

Tap End Session to complete programming session return to title screen

4

SCREEN INITIAL FINAL

STIMULATION UNDO **2** **3** GO TO STIMULATION

STIMULATION On On

GROUP A

Left ANT

2.3 V
90 μs
130 Hz

Right ANT

2.7 V
60 μs
130 Hz

Left ANT

2.3 V
90 μs
130 Hz

Right ANT

2.7 V
60 μs
130 Hz

CLINICIAN NOTES

NO ALERTS CLOSE END SESSION

Stimulation On John Smith 3:17 PM

INSIGHTFUL
INFORMATION
YOU NEED

The image displays two overlapping screenshots of the Medtronic mobile application. The top screenshot shows a menu with the following items: Usage, Reports, About, and Preferences. The bottom screenshot shows a 'SELECT TASK' screen with a list of tasks: SETUP, STIMULATION, IMPEDANCE, REPLACEMENT, and END SESSION. The interface includes a 'Stimulation' toggle switch set to 'On' and a 'NO ALERTS' button at the bottom left.

MENU OPTIONS

PATIENT
Patient Name: John Smith
Patient ID: MRN:123456
Diagnosis: Epilepsy

DEVICE
Device Model: Activa PC
Model Number: 37601
Serial Number: NKM_FOLLOWUP_PC
Implant Date: Apr 22, 2017
Battery Level: 3.81 V | OK
Cycling: Enabled
On Duration: 1 Minute
Off Duration: 5 Minutes

IMPEDANCE
Status: INVESTIGATE
LEFT HEMISPHERE

NO ALERTS

SELECT TASK

SETUP

Medtronic

- Usage
- Reports**
- About
- Preferences

SELECT TASK

- SETUP
- STIMULATION
- IMPEDANCE
- REPLACEMENT
- END SESSION

INSIGHTFUL INFORMATION YOU NEED - REPORTS

Various report
types


Session Report
[No name]
Activa PC NKM728744

Medtronic
Session Date:
Oct 24, 2018 1:35 PM

[No name]

Diagnosis:
Patient ID:

DEVICE

	NEUROSTIMULATOR	Serial Number	Implanted
	Activa PC	NKM728744	Oct 18, 2018
	Model Number	Firmware Version	
	37601	2.20	
	Neurostimulator Location	Battery Level	
	Chest Left	2.95 V OK	

STIMULATION SETTINGS


	Initial		Final	
	Stimulation: Off		Stimulation: Off	
	Left STN	Right STN	Left STN	Right STN
A	0.0 V 60 µs 125 Hz	0.0 V 60 µs 125 Hz	0.0 V 60 µs 125 Hz	0.0 V 60 µs 125 Hz
B	0.0 V 60 µs 125 Hz	0.0 V 60 µs 125 Hz	0.0 V 60 µs 125 Hz	0.0 V 60 µs 125 Hz
C	0.0 V 60 µs 130 Hz		0.0 V 60 µs 130 Hz	
D	0.0 V 60 µs 130 Hz		0.0 V 60 µs 130 Hz	

Session Report
[No name]
Activa PC NKM728744

Medtronic
Session Date:
Oct 24, 2018 1:35 PM

Clinician Notes

~:0/5 (0)
~X:5/5 (2)




A1	Session Date	A	B	C	D	E
	Session Date	Patient Name	Patient ID	Diagnosis	INS Model	INS SN
	Tue Oct 10 15:12:31 CDT 2017	John Smith	MRN:123456	Parkinson's Disease	37612	NKG_FOLLO

GOAL OF DBS PROGRAMMING


Deliver the therapy to the brain target of interest while minimizing stimulation of surrounding structures

- Using the lead electrode closest to the desired target provides maximal benefit and minimizes stimulation-induced adverse effects
- Setting appropriate stimulation parameters ensures that the desired brain target, but not adjacent structures, receives the stimulation



DBS PROGRAMMING ALGORITHM SAMPLE

INITIAL PROGRAMMING



FOLLOW-UP PROGRAMMING VISITS

- **Review interim changes**
 - Seizure status (review event counter)
 - Medication changes
 - Adverse effects
 - Medication induced
 - Stimulation induced
- **Interrogate device**
- **Assess where the stimulation parameters are within the therapeutic range**
- **Formulate and implement management plan**
 - Adjust stimulation
 - Adjust medication
 - Make no adjustments
- **End session , then synchronize patient programmer to reset their seizure counter**

INS SETTINGS (HEAD-ONLY AND FULL-BODY SCAN)

	type	Settings
Full Body Eligible	37602, 37603, 37601 (no pocket adaptor implanted)	Unipolar configuration ^a — Turn therapy off. Bipolar configuration ^b — Keep therapy on or turn therapy off.
	37603, 37601 (pocket adaptor implanted)	Turn therapy off.
Head Scan ONLY	37601	Turn therapy off.
	7428 KINETRA	Turn therapy off. Disable magnetic (reed) switch. Disable day cycling.
	7426 SOLETRA	Turn therapy off. Set to bipolar configuration. Set amplitude to 0 volts.

PC with NO pocket Adaptor

Full Body Eligible

PC WITH a pocket Adaptor

Head Scan ONLY

Bipolar/Unipolar Defined

^a At least one electrode is negative and the case is positive, the other electrodes can be either negative or off.
^b At least one electrode is positive, one electrode is negative, and the case is off.

MRI ELIGIBILITY SHEET (HEAD-ONLY/FULL-BODY)

EN MRI Eligibility Sheet

Completion of these sections is required prior to conducting an MRI scan on patients with a Medtronic Deep Brain Neurostimulator or implanted leads.

This form is for Medtronic DBS Systems only. Refer to the *MRI guidelines for Medtronic Deep Brain Stimulation Systems instructions for use* for complete safety information and instructions for conducting an MRI scan. The manual and the MRI eligibility sheet are available through www.medtronic.com/mri or by contacting Medtronic.

Complete or review the following sections to identify MRI scan-type eligibility (page 1) and to ensure that the DBS system is prepared for the MRI scan (page 2).

1. **Date eligibility was determined**

Patient name _____ Date of birth _____

DBS managing clinician _____ Clinician phone number _____

The MRI scan-type eligibility is based on the following combination of factors pertaining to the patient's DBS system:

- implantable neurostimulator model number
- presence of an implanted pocket adaptor
- lead implant status
- neurostimulation system integrity

2. **Medtronic DBS System types**

Implanted neurostimulation system

Lead-only system

3. This section is not applicable because the patient does not have an implanted neurostimulation system.

Medtronic DBS Neurostimulators	
Head-only eligible	Full-body eligible
<input type="checkbox"/> Model 37602 (Activa® SC)	<input type="checkbox"/> Model 37612 (Activa® RC)
<input type="checkbox"/> Model 7428 (Kinetra®)	<input type="checkbox"/> Model 37603 (Activa® SC)
<input type="checkbox"/> Model 7426 (Soletra®)	<input type="checkbox"/> Model 37601 (Activa® PC)

4. This section is not applicable because the patient does not have an implanted neurostimulation system.

Pocket adaptors	
Head-only eligible	Full-body eligible
<input type="checkbox"/> Pocket adaptor(s) implanted (with Activa models 37612, 37603, and 37601 neurostimulators)	<input type="checkbox"/> No pocket adaptor
<input type="checkbox"/> No pocket adaptor (with Activa SC Model 37602, Soletra Model 7426 or Kinetra Model 7428 neurostimulator)	

5. This section is not applicable because the patient does not have a lead-only system.

Lead-only systems	
Head-only eligible	Full-body eligible
<input type="checkbox"/> Partially-implanted lead(s)	<input type="checkbox"/> Fully-implanted lead(s)

6. Indicate the MRI scan-type eligibility based on the information provided. If more than one DBS system is implanted, MRI scan eligibility should be based on the most restrictive DBS system components.

Head-only MRI scan Full-body MRI scan

English

7. Verified system integrity for all implanted neurostimulators by testing for open and short circuits using the clinician programmer to measure impedance.

- This section is not applicable because the patient does not have an implanted neurostimulation system.
- System integrity verified. (No open or short circuits detected.)
- System is compromised. (Open or short circuit verified.) **DO NOT** perform the MRI.

8. Neurostimulator programmed to the recommended settings listed below. Managing clinicians should document settings before changing to help ensure accurate re-programming after the MRI scan. Ensure that section 3 and section 4 have been completed on page 1 of this eligibility sheet. If the sections have not been completed, an MRI scan should not be performed.

- This section is not applicable because the patient does not have an implanted neurostimulation system.

Recommended neurostimulator settings for MRI	
<input type="checkbox"/> 37612, 37603, 37601 (no pocket adaptor implanted)	<input type="checkbox"/> Therapy off in unipolar or bipolar configuration. <input type="checkbox"/> Therapy on in bipolar configuration.
<input type="checkbox"/> 37612, 37603, 37601 (with pocket adaptor implanted)	Therapy off.
<input type="checkbox"/> 37602	Therapy off.
<input type="checkbox"/> 7428	Therapy off. Magnetic (reed) switch disabled. Day cycling disabled.
<input type="checkbox"/> 7426	Therapy off. Bipolar configuration setting. Amplitude set to 0 volts.

- Neurostimulator has been programmed by clinician.

Radiology staff: Ask the patient if therapy settings have changed prior to the MRI appointment. If the patient has changed therapy settings, contact patient's clinician or Medtronic Technical Services.

- Neurostimulator will be programmed at MRI appointment by the patient or the clinician.

Radiology staff: Have the patient or clinician program the neurostimulator(s) now according to the recommended settings.

9. For lead-only systems, partially-implanted or fully-implanted leads have been prepared by the clinician for the MRI scan. Ensure that section 5 has been completed on page 1 of the eligibility sheet. If that section has not been completed, an MRI scan should not be performed.

- This section is not applicable because the patient does not have a lead-only system.
- Fully-implanted lead(s) has been capped and internalized.
- Partially-implanted lead(s) has been insulated, and the externalized portion of the lead is out-of-contact with the patient, is straight with no loops, and is centered in the head coil.

10. DBS clinician signature _____ Date _____

11. **FOR RADIOLOGY STAFF ONLY**

This form reviewed by _____

Print name _____ Signature _____

MRI Technologist/Radiographer Radiologist

Date _____

Assess other implanted medical devices before conducting an MRI scan.

SUMMARY

POST-OPERATIVE MRI IN PATIENTS WITH BOTH DBS & VNS

Post-op MRI with both DBS and VNS implants *may* be possible provided all safety conditions are met. Several questions must be considered prior to scanning:

1. **Determine what model of DBS system** is implanted and if full-body or head-only MRI eligible
2. **Determine what brand/model of VNS system** is implanted
3. Determine if the VNS is complete, active or inactive, or partially excised
4. If the VNS system is complete (leads & IPG), determine if VNS IPG is to be excised or left in place prior to MRI scan
5. If the **VNS IPG has previously been removed**, determine:
 - a. **If > 2 cm of VNS lead length remains, or**
 - b. **If < 2 cm of VNS lead length remains**

For MRI compatibility questions about any implanted VNS system or partial system, together with a DBS system, Review **both sets of MRI instructions**, and to follow the more restrictive instructions.

Having answers to the questions above will facilitate this process.

DISCUSSION QUESTIONS

- **How can we predict responders?**
- **How can we best target the implant?**
- **Why does benefit increase over time?**
- **Which parameters are best?**
- **Where does DBS fit in with other therapies?**

Discussion – Cases

BRIEF STATEMENT: MEDTRONIC DBS THERAPY FOR PARKINSON'S DISEASE, TREMOR, DYSTONIA AND EPILEPSY

Medtronic DBS Therapy for Parkinson's Disease, Tremor, Dystonia, and Epilepsy: Product labeling must be reviewed prior to use for detailed disclosure of risks.

Indications:

Medtronic DBS Therapy for Parkinson's Disease: Bilateral stimulation of the internal globus pallidus (GPi) or the subthalamic nucleus (STN) using Medtronic DBS Therapy for Parkinson's Disease is indicated for adjunctive therapy in reducing some of the symptoms in individuals with levodopa-responsive Parkinson's disease of at least 4 years' duration that are not adequately controlled with medication, including motor complications of recent onset (from 4 months to 3 years) or motor complications of longer-standing duration.

Medtronic DBS Therapy for Tremor: Unilateral thalamic stimulation of the ventral intermediate nucleus (VIM) using Medtronic DBS Therapy for Tremor is indicated for the suppression of tremor in the upper extremity. The system is intended for use in patients who are diagnosed with essential tremor or parkinsonian tremor not adequately controlled by medications and where the tremor constitutes a significant functional disability.

Medtronic DBS Therapy for Dystonia: Unilateral or bilateral stimulation of the internal globus pallidus (GPi) or the subthalamic nucleus (STN) using Medtronic DBS Therapy for Dystonia is indicated as an aid in the management of chronic, intractable (drug refractory) primary dystonia, including generalized and/or segmental dystonia, hemidystonia, and cervical dystonia (torticollis), in patients seven years of age or above.

Medtronic DBS Therapy for Epilepsy: Bilateral stimulation of the anterior nucleus of the thalamus (ANT) using the Medtronic DBS System for Epilepsy is indicated as an adjunctive therapy for reducing the frequency of seizures in individuals 18 years of age or older diagnosed with epilepsy characterized by partial-onset seizures, with or without secondary generalization, that are refractory to three or more antiepileptic medications.

The Medtronic DBS System for Epilepsy has demonstrated safety and effectiveness for patients who average six or more seizures per month over the three most recent months prior to implant of the DBS system (with no more than 30 days between seizures). The Medtronic DBS System for Epilepsy has not been evaluated in patients with less frequent seizures.

Contraindications: Medtronic DBS Therapy is contraindicated for patients who are unable to properly operate the neurostimulator and, for Parkinson's disease and Essential Tremor, patients for whom test stimulation is unsuccessful. The following procedures are contraindicated for patients with DBS systems: diathermy (e.g., shortwave diathermy, microwave diathermy or therapeutic ultrasound diathermy), which can cause neurostimulation system or tissue damage and can result in severe injury or death; Transcranial Magnetic Stimulation (TMS); and certain MRI procedures using a full body transmit radio-frequency (RF) coil, a receive-only head coil, or a head transmit coil that extends over the chest area if they have an implanted Soletra™ Model 7426 Neurostimulator, Kinetra™ Model 7428 Neurostimulator, Activa™ SC Model 37602 Neurostimulator, or Model 64001 or 64002 pocket adaptor.

BRIEF STATEMENT: MEDTRONIC DBS THERAPY FOR PARKINSON'S DISEASE, TREMOR, DYSTONIA AND EPILEPSY

Warnings and Precautions: There is a potential risk of brain tissue damage using stimulation parameter settings of high amplitudes and wide pulse widths and, for Parkinson's disease and essential tremor, a potential risk to drive tremor using low frequency settings. Extreme care should be used with lead implantation in patients with an increased risk of intracranial hemorrhage. Sources of electromagnetic interference (EMI) may cause device damage or patient injury. Theft detectors and security screening devices may cause stimulation to switch ON or OFF and may cause some patients to experience a momentary increase in perceived stimulation. The DBS System may be affected by or adversely affect medical equipment such as cardiac pacemakers or therapies, cardioverter/ defibrillators, external defibrillators, ultrasonic equipment, electrocautery, or radiation therapy. MRI conditions that may cause excessive heating at the lead electrodes which can result in serious injury, including coma, paralysis, or death, or that may cause device damage, include: neurostimulator implant location other than pectoral and abdominal regions; unapproved MRI parameters; partial system explants ("abandoned systems"); misidentification of neurostimulator model numbers; and broken conductor wires (in the lead, extension or pocket adaptor). The safety of electroconvulsive therapy (ECT) in patients receiving DBS Therapy has not been established. Tunneling the extension too superficially or too deeply may result in nerve or vascular injury, or tunneling through unintended anatomy. The lead-extension connector should not be placed in the soft tissues of the neck due to an increased incidence of lead fracture. Abrupt cessation of stimulation should be avoided as it may cause a return of disease symptoms, in some cases with intensity greater than was experienced prior to system implant ("rebound" effect). Onset of status dystonicus, which may be life-threatening, may occur in dystonia patients during ongoing or loss of DBS therapy. For Epilepsy, cessation, reduction, or initiation of stimulation may potentially lead to an increase in seizure frequency, severity, and new types of seizures. For Epilepsy, symptoms may return with an intensity greater than was experienced prior to system implant, including the potential for status epilepticus. Patients using a rechargeable neurostimulator for Parkinson's disease or Essential Tremor should check for skin irritation or redness near the neurostimulator during or after recharging, and contact their physician if symptoms persist. Loss of coordination in activities such as swimming may occur. Depression, suicidal ideations and suicide have been reported in patients receiving Medtronic DBS Therapy for Movement Disorders and Epilepsy, although no direct cause-and-effect relationship has been established. For Epilepsy, preoperatively, assess patients for depression and carefully balance this risk with the potential clinical benefit. Postoperatively, monitor patients closely for new or changing symptoms of depression and manage these systems appropriately. Patients should be monitored for memory impairment. Memory impairment has been reported in patients receiving Medtronic DBS Therapy for Epilepsy, although no direct cause and effect relationship has been established. The consequences of failing to monitor patients are unknown. When stimulation is adjusted, monitor patients for new or increased seizures, tingling sensation, change in mood, or confusion.

Adverse Events: Adverse events related to the therapy, device, or procedure can include intracranial hemorrhage, cerebral infarction, CSF leak, pneumocephalus, seizures, surgical site complications (including pain, infection, dehiscence, erosion, seroma, and hematoma), meningitis, encephalitis, brain abscess, cerebral edema, aseptic cyst formation, device complications (including lead fracture and device migration) that may require revision or explant, extension fibrosis (tightening or bowstringing), new or exacerbation of neurological symptoms (including vision disorders, speech and swallowing disorders, motor coordination and balance disorders, sensory disturbances, cognitive impairment, and sleep disorders), psychiatric and behavioral disorders (including psychosis and abnormal thinking), cough, shocking or jolting sensation, ineffective therapy, and weight gain or loss.

Safety and effectiveness has not been established for patients with previous surgical ablation procedures, dementia, coagulopathies, or moderate to severe depression, or patients who are pregnant. Parkinson's disease and essential tremor: safety and effectiveness has not been established for patients under 18 years or patients with neurological disease other than idiopathic Parkinson's disease or Essential Tremor. Essential tremor: safety and effectiveness has not been established for bilateral stimulation or for patients over 80 years of age. Dystonia: age of implant is suggested to be that at which brain growth is approximately 90% complete or above. Epilepsy: the safety and effectiveness of this therapy has not been established for patients without partial-onset seizures, patients who are pregnant or nursing, patients under the age of 18 years, patients with coagulopathies, and patients older than 65 years.

Humanitarian Device (Dystonia): Authorized by Federal Law as an aid in the management of chronic, intractable (drug refractory) primary dystonia, including generalized and/or segmental dystonia, hemidystonia, and cervical dystonia (torticollis), in patients seven years of age or above. The effectiveness of the devices for treating these conditions has not been demonstrated.

USA Rx only Rev 08/18

THANK YOU