


3T STIR (Short T1 Inversion Recovery)

Left: Unlabeled 3T STIR Sagittal Image

Right: Labeled 3T STIR Sagittal Image

(mtt = mammillothalamic tract, Apr + AM = ANT)

Source: Lehtimäki K, Mottonen T, Jarventausta K, Katisko J, Tahtinen T, Haapasalo J, Niskakangas T, Kiekara T, Ohman J, Peltola J. Outcome based definition of the anterior thalamic deep brain stimulation target in refractory epilepsy. *Brain stim.* 2016; 9: 268-275




1.5T STIR (Short T1 Inversion Recovery)

Labeled 1.5T STIR Coronal Image

(mtt = mammillothalamic tract, ANT = Anterior nucleus of the thalamus, eml = external medullary lamina)

Source: Jiltsova E, Mottonen T, Fahlstrom M, Haapasalo J, Tahtinen T, Peltola J, Ohman J, Larsson EM, Kiekara T, Lehtimäki K. Imaging of anterior nucleus of the thalamus using 1.5 T MRI for deep brain stimulation targeting in refractory epilepsy. *Neuromodulation.* 2016; 19(8):812-817



3T MPRAGE (Magnetization-Prepared Rapid Gradient Echo)

Sagittal Image:

Arrow indicates location of ANT

Arrowhead points to mammillothalamic tract

Source: Buentjen L, Kopitzki K, Schmitt FC, Voges J, Tempelmann C, Kaufmann J, Kanowski M. Direct targeting of the thalamic anteroventral nucleus for deep brain stimulation by T1-weighted magnetic resonance imaging at 3T. *Stereotact Funct Neurosurg.* 2014; 92:25-30



3T FGATIR (Fast Gray Matter Acquisition T1 Inversion Recovery)

Sagittal Image:

Red targeting dot placed within ANT

Arrow points to mammillothalamic tract

Source: Medtronic



Brief Statement: Medtronic DBS Therapy for Epilepsy

Medtronic DBS systems are MR Conditional which means they are marked to indicate they are safe in the MR environment as long as certain conditions are met. Read and fully understand the *MRI Guidelines for Medtronic deep brain stimulation systems* before conducting the MRI examination. Go to www.medtronic.com/mri or contact Medtronic at 1-800-707-0933 for a copy. Also review current MRI manufacturer labeling before conducting the MRI.

Medtronic DBS Therapy for Epilepsy: Product labeling must be reviewed prior to use for detailed disclosure of risks.

Indications: Bilateral stimulation of the anterior nucleus of the thalamus (ANT) using the Medtronic DBS System for Epilepsy is indicated as an adjunctive therapy for reducing the frequency of seizures in individuals 18 years of age or older diagnosed with epilepsy characterized by partial-onset seizures, with or without secondary generalization, that are refractory to three or more antiepileptic medications.

The Medtronic DBS System for Epilepsy has demonstrated safety and effectiveness for patients who average six or more seizures per month over the three most recent months prior to implant of the DBS system (with no more than 30 days between seizures). The Medtronic DBS System for Epilepsy has not been evaluated in patients with less frequent seizures.

Contraindications: Medtronic DBS Therapy is contraindicated for patients who are unable to properly operate the neurostimulator. The following procedures are contraindicated for patients with DBS systems: diathermy (e.g., shortwave diathermy, microwave diathermy or therapeutic ultrasound diathermy), which can cause neurostimulation system or tissue damage and can result in severe injury or death; Transcranial Magnetic Stimulation (TMS); and certain MRI procedures using a full body transmit radio-frequency (RF) coil, a receive-only head coil, or a head transmit coil that extends over the chest area if the patient has an implanted Soletra™ Model 7426 Neurostimulator, Kinetra™ Model 7428 Neurostimulator, Activa™ SC Model 37602 Neurostimulator, or Model 64001 or 64002 pocket adaptor.

Warnings and Precautions: There is a potential risk of brain tissue damage using stimulation parameter settings of high amplitudes and wide pulse widths. Extreme care should be used with lead implantation in patients with a heightened risk of intracranial hemorrhage. Sources of electromagnetic interference (EMI) may cause device damage or patient injury. Theft detectors and security screening devices may cause stimulation to switch ON or OFF and may cause some patients to experience a momentary increase in perceived stimulation. The DBS System may be affected by or adversely affect medical equipment such as cardiac pacemakers or therapies, cardioverter/defibrillators, external defibrillators, ultrasonic equipment, electrocautery, or radiation therapy. MRI conditions that may cause excessive heating at the lead electrodes which can result in serious injury, including coma, paralysis, or death, or that may cause device damage, include: neurostimulator implant location other than pectoral and abdominal regions; unapproved MRI parameters; partial system explants ("abandoned systems"); misidentification of neurostimulator model numbers; and broken conductor wires (in the lead, extension or pocket adaptor). The safety of electroconvulsive therapy (ECT) in patients receiving DBS Therapy has not been established. Tunneling the extension too superficially or too deeply may result in nerve or vascular injury, or tunneling through unintended anatomy. The lead-extension connector should not be placed in the soft tissues of the neck due to an increased incidence of lead fracture. Cessation, reduction, or initiation of stimulation may potentially lead to an increase in seizure frequency, severity, and new types of seizures. Symptoms may return with an intensity greater than was experienced prior to system implant, including the potential for status epilepticus. Loss of coordination in activities such as swimming may occur. Depression, suicidal ideations and suicide have been reported in patients receiving Medtronic DBS Therapy for Epilepsy, although no direct cause-and-effect relationship has been established. Preoperatively, assess patients for depression and carefully balance this risk with the potential clinical benefit. Postoperatively, monitor patients closely for new or changing symptoms of depression and manage these symptoms appropriately. Patients should be monitored for memory impairment. Memory impairment has been reported in patients receiving Medtronic DBS Therapy for Epilepsy, although no direct cause and effect relationship has been established. The consequences of failing to monitor patients are unknown. When stimulation is adjusted, monitor patients for new or increased seizures, tingling sensation, change in mood, or confusion.

Adverse Events:

Adverse events related to the therapy, device, or procedure can include intracranial hemorrhage, cerebral infarction, CSF leak, pneumocephalus, seizures, surgical site complications (including pain, infection, dehiscence, erosion, seroma, and hematoma), meningitis, encephalitis, brain abscess, cerebral edema, aseptic cyst formation, device complications (including lead fracture and device migration) that may require revision or explant, extension fibrosis (tightening or bowstringing), new or exacerbation of neurological symptoms (including vision disorders, speech and swallowing disorders, motor coordination and balance disorders, sensory disturbances, cognitive impairment, and sleep disorders), psychiatric and behavioral disorders (including psychosis and abnormal thinking), cough, shocking or jolting sensation, ineffective therapy and weight gain or loss.

The safety and effectiveness of this therapy has not been established for patients without partial-onset seizures, patients who are pregnant or nursing, patients under the age of 18 years, patients with coagulopathies, and patients older than 65 years.

USA Rx Only Rev 06/18